



L.H. SECTIONAL ELEVATION




FRONT ELEVATION


Recommendation- Socket to feed fan unit by others, dedicated circuit with suitable RCD protection for motors.




R.H. SECTIONAL ELEVATION



PLAN VIEW



SECTIONAL PLAN VIEW



FUME CUPBOARD CONTROL PANEL FASCIA

Exterior	16swg zinc coated mild steel, epoxy powder coated in textured white RAL 9003.
Interior	6mm thick clear toughened glass sides and 4.5mm thick white PVC back panel.
Baffles	4.5mm white PVC baffle.
Cill	16swg zinc coated mild steel, epoxy powder coated in textured white RAL9003
Worktop	Grey GRP complete with raised edges.
Support	38mm sq steel, fully welded with adjustable feet. Epoxy powder coated in textured white RAL9003.
Storage	Mild steel underbench unit with lockable doors, epoxy powder coated textured RAL9003 white.
Drainage	Vulcathene 501 Dripcup c/w 1.5" BSP male outlet. Termination of drainage by OTHERS
Sash	6mm clear toughened glass with stop @400mm that automatically resets
Light	LED light panel, over a sealed toughened glass lens
Services	FAR. Pre-plumbed to low level, for connection to mains by OTHERS
Electrics	Pre-wired to junction box, mounted on soffit. Supplied with 2 No. flying leads and plugs for connection to sockets provided on wall by OTHERS.
Controls	Electronic Control Panel (See above)

SEE SHEET 2 FOR OTHER SIZE OPTIONS

Holliday Fielding Hocking Limited
 Unit 9A, Aire Place Mills
 143 Kirkstall Road
 Leeds
 LS31JL
 Tel: 01132454111
 Fax: 0113 246 9964
 EMail: info@labfumecupboards.co.uk
 Web: www.labfumecupboards.co.uk

B	07/09/2018	PW	BROEN PARTS REPLACED BY FAR	IH
REV.	DATE	BY	DESCRIPTION	CHK

CLIENT:

DRAWING TITLE:
 EDUCATION FIXED DUCTED FUME CUPBOARD 1000(W) (THREE PHASE) SALES DRAWING

PART / ASSEMBLY NO.: P-0033-ASM-800 PROJECTION:


MATERIAL: SEE INSET

FINISH: SEE INSET


DRAWN BY:	SCALE:	PAPER:	ORIG. DATE:
PW	1:15	A3	31/05/2017

CHECKED:	SHEET:	PROJECT REF:
IH	1 of 2	


DRAWING NO.:	REVISION:
P-0033-ASM-819	B




1000







1200



1500



1800



Holliday Fielding Hocking Limited

Unit 9A, Aire Place Mills
143 Kirkstall Road
Leeds
LS31JL

Tel: 01132454111

Fax: 0113 246 9964

E-Mail: info@labfumecupboards.co.uk

Web: www.labfumecupboards.co.uk


DRAWING TITLE:

EDUCATION FIXED DUCTED FUME CUPBOARD
1000(W) (THREE PHASE) SALES DRAWING

PART / ASSEMBLY NO.:

P-0033-ASM-800

PROJECTION:



MATERIAL:

SEE INSET

FINISH:

SEE INSET

DRAWN BY: SCALE: PAPER: ORIG. DATE:

PW 1:20 A3 31/05/2017

CHECKED: SHEET: PROJECT REF:

IH 2 of 2

DRAWING NO.:

P-0033-ASM-819

REVISION:

B