



| | |
|-----------|---|
| Exterior | 14 SWG aluminium, powder coated textured white RAL 9003. |
| Interior | 6 mm thick clear toughened glass sides and back panel. |
| Baffles | 4 mm clear polycarbonate. |
| Cill | 14 SWG aluminium, powder coated textured white RAL 9003. |
| Worktop | Grey GRP complete with raised edges. |
| Storage | 1 x cupboard. |
| Drainage | Vulcathene 501 drip cup c/w 1.5" BSP male outlet. Flexible hose supplied for connection to drain. |
| Sash | 6 mm clear toughened glass. Counterbalanced with sash stop @360 mm that automatically resets. |
| Light | LED light panel, over a sealed toughened glass lens. |
| Services | FAR. Pre-plumbed to low level, for connection to docking station. |
| Electrics | Pre-wired to junction box, mounted on soffit. 240V 32A ring main connection to docking station. |
| Controls | Electronic Control Panel. |

Holliday Fielding Hocking Limited
 Unit 9A, Aire Place Mills
 143 Kirkstall Road
 Leeds
 LS31JL
 Tel: 01132454111
 Fax: 0113 246 9964
 EMail: info@labfumecupboards.co.uk
 Web: www.labfumecupboards.co.uk

| | | | |
|----------------------|--|--------------|-------------|
| CLIENT: | | | |
| DRAWING TITLE: | SSR2000 FUME CUPBOARD 1000(W) SALES DRAWING | | |
| PART / ASSEMBLY NO.: | P-0032-ASM-001 | PROJECTION: | |
| MATERIAL: | SEE INSET | | |
| FINISH: | SEE INSET | | |
| DRAWN BY: | SCALE: | PAPER: | ORIG. DATE: |
| PW | 1:15 | A3 | 12/07/2017 |
| CHECKED: | SHEET: | PROJECT REF: | |
| IH | 1 of 2 | | |
| DRAWING NO.: | P-0032-ASM-022 | | REVISION: |
| | | | E |



VIEW ON FRONT



VIEW ON REAR

POWER LEAD
c/w 13A PLUG

RESTRAINING CABLE

HANDLE

GAS, WATER AND WASTE OUTLET

Holliday Fielding Hocking Limited

Unit 9A, Aire Place Mills
143 Kirkstall Road
Leeds
LS31JL
Tel: 01132454111
Fax: 0113 246 9964
EMail: info@labfumecupboards.co.uk
Web: www.labfumecupboards.co.uk

CLIENT:

DRAWING TITLE:
SSR2000 FUME CUPBOARD 1000(W)
SALES DRAWING

| | |
|----------------------|-------------|
| PART / ASSEMBLY NO.: | PROJECTION: |
| P-0032-ASM-001 | |

MATERIAL:
SEE INSET

FINISH:
SEE INSET

| | | | |
|-----------|--------|--------|-------------|
| DRAWN BY: | SCALE: | PAPER: | ORIG. DATE: |
| PW | 1:10 | A3 | 12/07/2017 |

| | | |
|----------|--------|--------------|
| CHECKED: | SHEET: | PROJECT REF: |
| IH | 2 of 2 | |

| | |
|----------------|-----------|
| DRAWING NO.: | REVISION: |
| P-0032-ASM-022 | E |